

Minimum technical requirements recommended for the implementation of **NemoStudio Family 2D Software**

Operating System	Minimum: Windows 7, Windows 8.x, Windows 10 (professional version); 32 bits. Recommended: Windows 10 Professional 64 bits.
Processor	Minimum: Intel Core i3, i5. Recommended: Intel Core i7 at 2.8Ghz.
RAM Memory	Minimum: 2 Gb. Recommended: 4 Gb.
Graphic Card	Minimum: GPU card with minimum 1 Gb of total VRAM dedicated. Recommended: GPU card with 2 Gb of total dedicated VRAM. (NVIDIA Card is recommended as it is the most tested).
Card Disc	Between 3 and 5 Gb for the software installation. Additional space may be needed, depending on the volume of cases.
Firewall configuration	Allow access to TCP/IP communication to executables NemoServer.exe and NemoStudio.exe
Network Speed (wired or wireless)	Network local installation. Minimum 1Gb.
Internet Access	Yes, mandatory for licensing control.
Mouse	A mouse with three buttons and a mousewheel is highly recommended.
For Apple Macbook, iMac	Bootcamp (native Windows) with the highest performance. Virtual Pc running Windows in a virtual environment. Parallels is recommended. It is advisable to configure the virtual environment with the minimum Windows requirements aforementioned.

Minimum technical requirements recommended for the implementation of **NemoStudio Family 3D Software**

Operating System	Minimum: Windows 7, Windows 8.x, Windows 10 (professional version); 32 bits. Recommended: Windows 10 Professional 64 bits.
Processor	Minimum: Intel Core i7 at 2.8Ghz.
RAM Memory	Minimum: 4 Gb. Recommended: 8 Gb.
Graphic Card	Minimum: GPU card with minimum 2 Gb of total dedicated VRAM (NVIDIA Card is recommended as it is the most tested).
Card Disc	5 Gb for the software installation. Additional space may be needed, depending on the volume of cases.
Firewall configuration	Allow access to TCP/IP communication to executables NemoServer.exe and NemoStudio.exe
Network Speed (wired or wireless)	Network local installation. Minimum 1Gb.
Internet Access	Yes, mandatory for licensing control (virtual id).
Mouse	A mouse with three buttons and a mousewheel is highly recommended.
For Apple Macbook, iMac	Bootcamp (native Windows) with the highest performance. Virtual Pc running Windows in a virtual environment. Parallels is recommended. It is advisable to configure the virtual environment with the minimum Windows requirements aforementioned.