



## Technical Specifications for NemoStudio

### Cloud and Desktop

|   |  |
|---|--|
| Operating System  | <b>Requirements:</b> Native Windows 10 Professional 64 bits or Windows 10 Home (with all Windows updates installed).   |
| Processor   | <b>Minimum Requirements:</b> Intel Core i7 at 2.8 Ghz.   |
| RAM Memory  | <b>Minimum Requirements:</b> 16 GB.  |
| Graphics Card   | NVidia graphics card. Minimum 2 GB. VRAM   |
| Hard Disc   | <b>15 GB</b> for the software installation. Additional space may be needed, depending on the volume of cases.  |
| Internet Access and Firewall configuration                    | <b>Mandatory:</b> Allow access to the Internet for getting the registered licenses of the customer. Allow access to the LAN to communicate with the NemoStudio application.  |
| Running in local area network (LAN)                           | Network cards wired, minimum speed 1 Gb./sec.  |
| Minimum WIFI network connection requirements (only for Cloud) | Internet speed minimum 100 Mbps symmetrical (Recommended at least 300 Mbps symmetrical), available from the computer itself. (Only for NemoStudio Cloud)   |
| Mouse   | <b>Highly recommended</b> external mouse with 2 buttons and the third being a wheel.   |
| For Apple Macbook, iMac                                       | <b>Bootcamp (native Windows)</b> for the highest performance. All previous requirements for Operating System and Hardware apply in this configuration too, such as: <b>Minimum Requirements:</b> Intel Core i7 at 2.8 Ghz. |

Table 1: Technical Specifications for NemoStudio Cloud and Desktop