

Technical Requirements for NemoStudio Cloud and Desktop

Operating System	Requirements: Windows 10 or posterior 64-bit (for Windows version); Internet browser updated to the latest versions (for web version)
Time of the system	Mandatory requirement: The Windows "Settings - Time & Language - Date & Time" option for " Set the time automatically " is enabled"
Processor	Minimum Requirements: Intel Core i7 at 2.8GHz or superior
RAM Memory	Minimum Requirements: 16 GB.
Graphics Card	NVidia graphics card. Minimum 4 GB. VRAM
Hard Disc	15 GB for the software installation. Additional space may be needed, depending on the volume of cases.
Interface	Minimum Requirements: Screen Resolution 1440 x 900
Internet Access and Firewall configuration	Mandatory: Allow access to the Internet for getting the registered licenses of the customer. Allow access to the LAN to communicate with the NemoStudio application.
Minimum WIFI network connection requirements (only for Cloud)	Internet speed minimum 100 Mbps symmetrical (Recommended at least 300 Mbps symmetrical), available from the computer itself. (Only for NemoStudio Cloud)
Mouse	Highly recommended external mouse with 2 buttons and the third being a wheel.
For Apple Macbook, iMac	Bootcamp (native Windows) for the highest performance. All previous requirements for Operating System and Hardware apply in this configuration too, such as: Minimum Requirements: Intel Core i7 at 2.8 Ghz.
Connection to the internet	Internet connection is required.

Table 1: Technical Requirements for NemoStudio Cloud and Desktop

Technical Requirements for NemoStudio Web

Processor	Minimum Requirements: Intel Core i5 at 2.8 Ghz.
RAM Memory	Minimum Requirements: 8 GB.
Graphics Card	NVidia graphics card. Minimum 1 GB. VRAM
Interface Requirement	Minimum Requirements: Screen Resolution 1440 x 900
Internet Access and Firewall configuration	Mandatory: Allow access to the Internet for getting the registered licenses of the customer. Allow access to the LAN to communicate with the NemoStudio application.
Minimum WIFI network connection requirements	Internet speed minimum 100 Mbps symmetrical (Recommended at least 300 Mbps symmetrical), available from the computer itself.
Mouse	Highly recommended external mouse with 2 buttons and the third being a wheel.
Navegator	Chrome or Edge
Processor	Minimum Requirements: Intel Core i7 at 2.8GHz or superior

Table 1: Technical Requirements for NemoStudio Web